


PERMANENT VS TEMPORARY PRISMS


INCREASED FIELD EXPANSION

Permanent prisms offer 76% more field expansion than temporary prisms!



INCREASED COLLISION DETECTION

Permanent Prisms offer more than double the potential collision detection over temporary prisms. Walk confidently through a grocery store or a mall.



QUALITY

Permanent Prisms are made from a top quality PMMA acrylic, giving better contrast, and more visual field. And they're glued into the lens using a resilient UV-cured adhesive. They simply won't peel off or fall out.

ASK ABOUT PERMANENT PRISMS

Peripheral Prisms are most helpful when you're walking or driving. But really, how much of the time are you walking or driving? Chances are, not too much. But when you are, you need something quick to give you awareness of your blind field.

With that in mind, we have many combinations of options to help fit both your lifestyle and your budget.

The Peli Lens™ is a trademark of Chadwick Optical Inc.



PELI LENS™

REAL FIELD EXPANSION FOR HEMIANOPSIA


PATIENT TAKE HOME GUIDE



How to Get the Most
From Your Peripheral Prisms

PLACEMENT

It is very important that you look between the prisms to get the benefits of the expanded visual field. Looking directly through the prism will cause double vision. This is undesirable.




When first worn, peripheral prisms create a certain level of visual confusion. You try to distinguish what information is in your blind field, and what information is in your seeing field. The goal of training is to minimize this confusion.

HOW TO REAPPLY

In the event that a prism falls off your lens, use the tracing template below to reapply your training lens.

1. Align prisms on lens to prisms on template.
2. Trace around the edge of your frame.
3. Use the trace of your frame to reapply prism.



Important: Do not remove prisms from spectacle lens for cleaning. For proper care of your temporary prisms, use the below procedure often. If the prisms become discolored or permanently detached, see your eye care professional for replacement.

- 1 Rinse spectacle lenses under a gentle stream of warm running water. If contaminants remain after rinsing, use a brush to clean grooves.
- 2 Pat or blot dry with a soft, lint-free cloth.
- 3 Follow any special cleaning procedures as instructed by your eye care specialist.

OUR MOST POPULAR PERMANENT PRISM OPTION



Chadwick's SLAM technology is a single removable lens that attaches to your glasses lenses via magnets.

This makes it super easy for you to quickly and deftly attach the lens when you need field expansion, then remove it when you don't.

TRAINING WITH REACH AND TOUCH

These exercises should be practiced frequently during the temporary prism trial period. Your goal before your follow-up appointment should be six hours of training, or about 20 minutes per day.

With a Buddy: While you fixate on the nose of your buddy, they move their hand into your blind field. You point to and touch their hand as you detect it through the prism.

DIY Method: You can do this training by yourself by reaching and touching objects as you detect them through the prisms.

TRAINING WHILE WALKING

Walk from uncluttered areas such as a hallway to progressively cluttered areas such as a living room filled with furniture as potential obstacles. Note your observations. When doing this exercise at home or in public places, you may want someone to assist you.

